

MEDICAL CARRIERS

USER GUIDE

By Appointment to Her Majesty the Queen Supplier of Tamper Evident Bags and Security Seals Versapak International Limited, Erith, Kent

OUTWARD

OUTWARD

CONTENTS

Carrier

Thermal Barriers

750g Freeze Boards

Maintain temperature up to











MEDIUM INSULATED CARRIER









LARGE **INSULATED CARRIER**











You need to seal the bag to be UN3373 compliant. Use Versapak's T2 Security Seals

t: 020 8333 5300

e: sales@versapak.co.uk

w: versapak.co.uk



1. Place contents inside



2. Place Thermal Barriers



3. Place Freeze Boards



4. Close and Seal



