

Large Insulated Medical Carrier - Pathology



Keeps contents cool for 20 hours

PYTB3

The Versapak Large Insulated Pathology Carrier is used by Pathology departments for storing and carrying specimens and samples within hospitals and medical centres at stable temperatures. It meets the requirements of UN3373 and Packing Instruction P650. Once filled, it is sealed using one of our tamper evident security seals. If the seal is broken, the recipient can see that the integrity of the contents has been compromised.



PREVENTS
INFECTION



UN3373 & P650
COMPLIANT



5 YEAR GUARANTEE



NHS APPROVED
PRODUCTION



INSULATION
TECHNOLOGY

Contact Us

UK & Ireland

e: info@versapak.co.uk
w: versapak.co.uk
t: 020 8333 5300

Rest of the World Email

e: info@versapakinternational.com
w: versapakinternational.com
t: +44 (0)20 8333 5300

CATEGORY	Insulated Medical Carriers	
PRODUCT NAME	Versapak Large Insulated Medical Carrier - Pathology	
ITEM NUMBER	PYTB3	
KEY FEATURES	 <ul style="list-style-type: none"> ▪ Built-in rigid impact protectors and internal padding ▪ Maintains the temperature of the contents for up to 20 hours ▪ Incorporates the Versapak tamper evident zip closure for security seal ▪ Printed with UN3373 logo as standard ▪ Tamper evident and traceable, when used with numbered or barcoded security seals ▪ "Easy change" flip/flap address window for outward and return addresses ▪ 2 x A5 windows for site location and bag content 	
DIMENSIONS (WxHxD)	External: 650 x 380 x 255mm Internal: 590 x 300 x 195mm Internal Capacity: 34.5 litres (Gross)	
SPECIFICATION	<ul style="list-style-type: none"> ▪ Material - High quality PVC coated nylon polyester ▪ Colours: Available in Blue / Green as standard ▪ Optional Extras: Shoulder Strap, Freeze Boards, Thermal Barriers, Dataloggers, Seals ▪ Customisation and overprinting options available subject to a minimum order quantity - contact us for details 	
SECURITY-T2 SECURITY SEALS Various styles and prints. You need to seal the carrier to be UN3373 compliant.	RELATED ITEMS - FREEZE BOARD Various freeze boards to maintain internal temperatures, contact us for details	RELATED ITEMS - THERMAL BARRIER Versapak provide the thermal barrier to separate transported items from freeze boards
		
  		
<p>Versapak International Limited is registered in England, registered number 03566179.</p> <p>Our registered address is The Versapak Centre, 4 Veridion Way, Erith, Kent, DA18 4AL.</p>		